


Drosophila




Drosophila Vials


FEATURES

- * Vials are precision moulded from premium grade polystyrene
- * Glass-like clarity for easy sample viewing
- * Disposable vials are a safe and cost-saving alternative to glass vials

The tray pack is great for quick, effortless filling



Light weight and easy for handling



VIAL DESIGN

- * Available in narrow (25 x 95mm) and wide (28.5 x 95mm)
- * Both narrow and wide vials are available in bulk and tray packaging
- * The tray pack 10 x 10 row configuration is great for quick, effortless filling




If they fall, they will not break like glass vials do

You can put your logo here



51-0500



51-0800



ORDERING INFORMATION

Cat #	Description	Size	Pack	Sterile	Unit
51-0500	Narrow	25 × 95mm	Bulk	No	25 Pieces/bag, 500 Pieces/case
51-0501	Narrow	25 × 95mm	Tray	No	100 Pieces/rack, 500 Pieces/case
51-0800	Wide	28.5 × 95mm	Bulk	No	25 Pieces/bag, 500 Pieces/case
51-0801	Wide	28.5 × 95mm	Tray	No	100 Pieces/rack, 500 Pieces/case

PRODUCT PARAMETER

Cat #	Description	Dimensions=mm		
		A	B	C
51-0500	Drosophila Vial	95	O.D. 22.6	O.D. 25
51-0800	Drosophila Vial	95	O.D. 25.4	O.D. 28.5

*See Biologix Drosophila Plugs on page 52.

Drosophila Bottles

FEATURES

- * Made of polypropylene and can be autoclaved, yet are economical enough to throw away after single use
- * A pliable top permits easy closure removal
- * The square-bottom is also transparent and features a firm, stiff neck for easier handling



BOTTLE DESIGN

- * 6oz (177ml) bottle
- * Molded graduations in 25ml increments up to 100ml
- * Square bottom
- * Available in both bulk and tray packaging




ORDERING INFORMATION

Cat #	Pack	Sterile	Unit
51-17720	Bulk Pack	No	200 Pieces/bag, 200 Pieces/case
51-17721	Rack Pack	No	25 Pieces/rack, 200 Pieces/case

PRODUCT PARAMETER

Cat #	Description	Dimensions=mm							
		A	B	C	D	E	F	G	H
51-17720	Drosophila Bottle	103.3	74.6	12.0	6.0	0.D.43.0	0.D.38.5	55.6	34.2

*See Biologix Drosophila Plugs on page 52.



Drosophila Plugs

FEATURES

- * Made of high density foam
- * Offer minimal shedding as well as reduce evaporation
- * Superior barrier to mite infestation
- * Exact fit with narrow 0.9 in. (25mm) or wide 1.1 in. (28.5mm) diameter vials or 6oz (177ml) bottles
- * Stretches yet maintains tight fit after repeated removal and closing
- * Flat surface accepts marking pen




Flat surface accepts marking pen

Fit well with bottles



Fit well with Vials



ORDERING INFORMATION

Cat #	Applied #	Fits	Sterile	Unit
51-0505	51-0500/0501	Narrow Plastic Vials	No	100 Pieces/bag, 10 Bags/case
51-0805	51-0800/0801	Wide Plastic Vials	No	100 Pieces/bag, 10 Bags/case
51-17725	51-17720/17721	6oz Plastic Bottles	No	100 Pieces/bag, 10 Bags/case

*See Biologix Drosophila Vials and Bottles from page 50 to 51.