

Gloves

FEATURES



Powder Free Latex Disposable Gloves

- * Powder-free, single chlorinated to provide good dry hand donning capability
- * Full textured surface enhances grip for efficient instrument handling
- * Good grip to handle dry instruments when moderate amount of handling is required
- * Good Thickness Distribution for tactile sensitivity
- * Suitable for all types of medical examination, home care and diagnostic and laboratory use
- * Biocompatibility Tested (ISO 10993)
- * Ambidextrous
- * Non-Sterile
- * Natural to Light Yellow Color



Powder Free Nitrile Disposable Gloves

- * Non-added components (sulphur and accelerator chemicals) to minimize exposure to Type IV hypersensitivity
- * Non-latex formulation eliminates exposure to Type 1 latex protein sensitivity
- * Nitrile material promotes better puncture and chemical/ chemo resistance properties
- * Low Dermatitis Potential Claim
- * Finger textured with improved working grip to enable efficient instrument handling in wet and dry condition
- * Beaded cuff design to don gloves efficiently and securely
- * Extra Length
- * Biocompatibility Tested (ISO 10993)
- * Ambidextrous
- * Non-Sterile
- * Dark Blue Color



ORDERING INFORMATION

Cat #	Description	Size	Color	Sterile	Unit
97-1111	Latex	XS	Nature to Light Yellow	No	100 Pieces/pack, 10 Packs/case
97-1112	Latex	S	Nature to Light Yellow	No	100 Pieces/pack, 10 Packs/case
97-1113	Latex	M	Nature to Light Yellow	No	100 Pieces/pack, 10 Packs/case
97-1114	Latex	L	Nature to Light Yellow	No	100 Pieces/pack, 10 Packs/case
97-1115	Latex	XL	Nature to Light Yellow	No	100 Pieces/pack, 10 Packs/case
97-6111	Nitrile	XS	Dark Blue	No	100 Pieces/pack, 10 Packs/case
97-6112	Nitrile	S	Dark Blue	No	100 Pieces/pack, 10 Packs/case
97-6113	Nitrile	M	Dark Blue	No	100 Pieces/pack, 10 Packs/case
97-6114	Nitrile	L	Dark Blue	No	100 Pieces/pack, 10 Packs/case
97-6115	Nitrile	XL	Dark Blue	No	100 Pieces/pack, 10 Packs/case

