

Histology



Slide Storage Boxes

FEATURES

- * Made from strong ABS material which is unaffected by humidity and is thoroughly insect proof
- * Store 1" x 3" (25mm x 75mm) microscope slides
- * Superior design prevents slides breakage and block dust from contaminating slides
- * Each slot is numbered to correspond to the slide inventory sheet on the inside cover
- * The 100-place box is secured with a rust-resistant clasp and hinge pin




BOX DESIGN

- * Available in 25-place or 100-place volume
- * Choose from either traditional cork or foam
- * 5 assorted colors (red, yellow, green, blue, and white) make the boxes easy to identify
- * Boxes are stackable to save space



- This is easy to be labeled
- Easy for each sample notes
- Rust-resistant clasp to secure the box
- Stackable to save space



ORDERING INFORMATION

Cat #	Description	Lining	Place	Sterile	Unit
41-5025	Assorted Colors	Foam	25-place	No	20 Pieces/case
41-5100	Assorted Colors	Foam	100-place	No	20 Pieces/case
41-6025	Assorted Colors	Cork	25-place	No	20 Pieces/case
41-6100	Assorted Colors	Cork	100-place	No	20 Pieces/case


PRODUCT PARAMETER

Cat #	Description	Dimensions=mm								
		A	B	C	D	E	F	G	H	I
41-6025	Box	140.0	133.0	62.5	88.0	71.4	80.0	23.0	16.4	6.6
41-6025C	Lid	140.0	137.0	82.8	88.0	13.7	16.1	-	-	-

*See Biologix Microscope Slides on page 46.





Microscope Slides



FEATURES

- * The slides are with clean glass surface
- * Ground edges prevent personal injury and chipping of glass caused by the cut edges
- * The glass slide end faces are polished
- * The slides without frosted area can be effectively used
- * The slides with frosted area are convenient for classification and storage

Polished end prevents personal injury



With frosted end and pencils can write on the area. Since the written characters are indelible in staining fluids or solvents, the glass slides of this type are convenient for classification and storage.

ORDERING INFORMATION


Cat #	Material	Type	Size	Sterile	Unit
7101	Clear Glass	Plain Slides	25.4 x 76.2mm	No	72 Pieces/pack, 50 Packs/case
7105	Clear Glass	Single Frosted End	25.4 x 76.2mm	No	72 Pieces/pack, 50 Packs/case

*See Biologix Slide Storage Boxes and Slide Mailers on page 45 and 47.

Slide Mailers

FEATURES

- * Made of sturdy cardboard to protect slides
- * Finger slots for easy removal of slides
- * Folding lids protect samples on slides
- * Storage and a safe transport method of microscope slides



DESIGN

- * Available in 1-place, 2-place or 20-place volume
- * 20-place mailer includes an inventory sheet and storage box
- * For 20-place mailers, they are packed 10 pieces in one inner durable pack to better protect products




- Sturdy cardboard to protect slides
- Packed by durable box to better protect mailers
- Finger slots for easy removal of slides
- The writing area is easy for notes

ORDERING INFORMATION

Cat #	Place	Sterile	Unit
41-4001	1-place	No	100 Pieces/pack, 8 Packs/case
41-4002	2-place	No	100 Pieces/pack, 8 Packs/case
41-4020	20-place	No	10 Pieces/pack, 4 Packs/case

PRODUCT PARAMETER


Cat #	Description	Dimensions=mm				
		A	B	C	D	E
41-4020	Slide Mailers	415.0	340.0	200.0	26.0	82.0



*See Biologix Microscope Slides on page 46.

41-4020

Embedding Cassettes/O-rings




FEATURES

- * Materials: high density polymer
- * Extreme chemical resistance



ORDERING INFORMATION

Cat #	Name	Color	Holes	Sterile	Unit
41-1001	Histology Cassettes	White	Strip	No	500 Pieces/pack, 2 Packs/case
41-1002	Histology Cassettes	Pink	Strip	No	500 Pieces/pack, 2 Packs/case
41-1003	Histology Cassettes	Green	Strip	No	500 Pieces/pack, 2 Packs/case
41-1004	Histology Cassettes	Yellow	Strip	No	500 Pieces/pack, 2 Packs/case
41-1005	Histology Cassettes	Blue	Strip	No	500 Pieces/pack, 2 Packs/case
41-2001	Biopsy Cassettes	White	Square	No	500 Pieces/pack, 2 Packs/case
41-2002	Biopsy Cassettes	Pink	Square	No	500 Pieces/pack, 2 Packs/case
41-2003	Biopsy Cassettes	Green	Square	No	500 Pieces/pack, 2 Packs/case
41-2004	Biopsy Cassettes	Yellow	Square	No	500 Pieces/pack, 2 Packs/case
41-2005	Biopsy Cassettes	Blue	Square	No	500 Pieces/pack, 2 Packs/case
41-3001	Embedding O-rings	White	-	No	500 Pieces/pack, 2 Packs/case
41-3002	Embedding O-rings	Pink	-	No	500 Pieces/pack, 2 Packs/case
41-3003	Embedding O-rings	Green	-	No	500 Pieces/pack, 2 Packs/case
41-3004	Embedding O-rings	Yellow	-	No	500 Pieces/pack, 2 Packs/case
41-3005	Embedding O-rings	Blue	-	No	500 Pieces/pack, 2 Packs/case



CASSETTE DESIGN

- * Cassettes available with strip or square holes
- * Feature a 45° angle writing area
- * Available in white, pink, green, yellow, and blue

PRODUCT PARAMETER

Cat #	Description	Dimensions=mm									
		A	B	C	D	E	F	G	H	I	J
41-1001	Cassette	28.8	27.0	35.5	31.1	27.0	29.0	81.6	30.4	45°	6.0